

Keep your soil tester accurate by using this simple cleaning procedure:

Place the Conditioning Film between the two metal plates of your soil pH tester (one on each side of the lower part of the housing). Squeeze lightly and rotate a few times.

This rubs the plates clean. They might have looked clean before, but they are often not chemically clean.

- Wipe plates with a clean rag or paper towel.
- Be sure not to touch plates with your fingers.
- Use the dull side of the film.
- One sheet is good for approximately 20 uses or until rubbing fails to cause the needle on the meter to deflect to the right when inserted in the soil.
- 10 sheets per pack.